

# clear55A<sup>®</sup>

methafilcon A lenses



Suddenly... everything is clear<sup>®</sup>

- + *Enhanced Comfort*
- + *Extended Wearing Time*
- + *Improved Ocular Health*
- + *Positive Wearer Preference*

clear55<sup>®</sup> was designed to provide an optimal balance of handling, comfort and vision. The lens has been enhanced with aberration correcting aspheric optics to optimize the overall visual experience and encompasses a design philosophy that addresses the wearer's handling and comfort needs.

-12.00  
+6.00

## OVERALL BALANCE



The performance of clear55A<sup>®</sup> has been extensively tested through a range of clinical studies including a 120 patient double masked crossover study where it was compared with other methafilcon lenses.

## OPTIMIZING VISUAL ACUITY THROUGH INNOVATIVE ASPHERIC DESIGN

clear55A<sup>®</sup> visual acuity is optimized over the power range and the incorporation of aberration correcting aspherics provides wearers with maximum contrast sensitivity, particularly in low illumination where it is needed most.

## OPTIMIZING WEARER COMFORT, FIT & HANDLING

clear55A<sup>®</sup> peripheral and edge design provides a balance between excellent comfort and assured handling. The comfort and handling of clear55A<sup>®</sup> was preferred by wearers compared with other leading brands, clear55A<sup>®</sup> larger diameter assists in centering the lens on the eye and ensures a latitude of fit over a wide range of corneal sizes, providing practitioners with a versatile fitting lens.



SPECIFICATIONS - clear55A <sup>®</sup>	
Material	45% methafilcon A, 55% water
Base Curve Diameter	8.7 mm 14.5 mm
Handling Tint	Light Blue
Power Range	-12.00 to +6.00 0.00 to -6.00 (0.25D steps) -6.50 to -12.00 (0.50D steps) +0.25 to +4.00 (0.25D steps) +4.50 to +6.00 (0.50D steps)
Centre Thickness*	0.09 mm @ -3.00D
Oxygen Permeability Dk @ 35°C	18.8 x 10 <sup>-11</sup> (cm <sup>2</sup> /s)(mlO <sub>2</sub> /ml.mm Hg)
Pack Size	6 lenses

\*-3.00 D lens

## POSITIVE WEARER PREFERENCE

The overall balance of clear55A<sup>®</sup> was appreciated by wearers who showed an overwhelming preference for the lens compared to other leading brands. The lens was preferred more often for its overall comfort and handling.

Convert your *low margin* contact lens practice into a *high margin* practice

O

**ocular insight**  
Suddenly... everything is clear<sup>®</sup>

phone: 877-ocular5 • fax: 877-ocular6 • email: info@ocularinsight.com • web: ocularinsight.com

Circle covers, month, date, phone, logo, name or graphic used or displayed in this brochure contain the registered and unregistered trademarks of Ocular Insight Inc. and Ocular Insight, Inc. All other trademarks indicated or registered on this brochure are registered in the territories in the United States and Singapore and may also be registered in other territories. Nothing on this brochure should be construed as granting, by implication or otherwise, any license or right to use, register, copy or reproduce any of the trademarks displayed in this brochure (in any form) without the prior written consent from Ocular Insight Inc. and Ocular Insight, Inc. All rights reserved.